










### 8.1 ΣΕΤ ΕΞΑΡΤΗΜΑΤΩΝ ΟΡΘΟΠΑΙΔΙΚΗΣ


| Περιγραφή  | Κωδικός        | Τεμάχια  |
|--|----------------|----------|
| <p><b>1.</b></p> <p><b>Pad for operation at intervertebral disc</b><br/>foam pad is electrically conductive;<br/>dimensions (LxWxH): 33 in x 18 in x 4 in<br/>(840 mm x 460 mm x 100 mm)</p>              | <b>1409494</b> | <b>2</b> |
| <p><b>2.α</b></p> <p><b>Adapter for head positioning</b><br/>to attach to TRUMPF coupling point system; for<br/>attaching head positioning accessories; stainless steel;<br/>spacer: 9 in (220 mm)</p>  | <b>4544785</b> | <b>2</b> |



| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <p><b>Skull clamp</b><br/>to attach to 1218969 skull clamp adapter; stable and precise three-point skull fixation; cast aluminum</p>    | 1218968 | 2       |
| <p><b>Basic unit neuro</b><br/>for attaching 1218969 Skull clamp adapter;<br/>simultaneous locking of longitudinal and vertical positions; used for prone and supine position of the patient; cast aluminum</p>  | 1218970 | 2       |

| Περιγραφή   | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Head rest for adults</b><br/>adjustable to head size; to connect to 1218969 skull clamp adapter; with detachable, elastic gel pad; used for prone and supine position of the patient;<br/>cast aluminum</p>      | 1218974 | 2       |
| <p><b>Head rest for children</b><br/>adjustable to head size; to connect to 1218969 skull clamp adapter; with detachable, elastic gel pad; used for prone and supine position of the patient;<br/>cast aluminum</p>  | 1218973 | 2       |



| Περιγραφή   | Κωδικός               | Τεμάχια         |
|---|-----------------------|-----------------|
| <p><b>Skull clamp adapter</b><br/>                     connection between 1218970 basic unit neuro, 1218968 skull clamp or 4145424 multi-purpose head rest; cast aluminum</p>  | <p><b>1218969</b></p> | <p><b>2</b></p> |
| <p><b>Skull pins for adults</b><br/>                     sterilisable; made of stainless steel; package with three pieces</p>   | <p><b>4145420</b></p> | <p><b>2</b></p> |


| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <b>Skull pins for children</b><br>sterilisable; made of stainless steel;<br>package with three pieces<br> | 4145421 | 2       |

|   |         |   |
|---|---------|---|
| <b>2β. Retainer for head rests</b><br>to attach to adapter for head positioning, shoulder<br>segment, shoulder adapter and ophthalmology adapter;<br>stainless steel<br> | 1317948 | 2 |
|---|---------|---|



| Περιγραφή   | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Connection holder basic asymmetric</b><br/>to attach to retainer for head rest; stainless steel</p>    | 1305449 | 2       |
| <p><b>Head rest basic</b><br/>with clamp; to attach to connection holder; Comfort Plus pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure; dimensions (LxV): 10 in x 8 in (250 mm x 200 mm)</p>  | 4549730 | 2       |




**3.**



| Περιγραφή   | Κωδικός               | Τεμάχια         |
|---|-----------------------|-----------------|
| <p><b>Supporting pad universal</b><br/>                     pad for extension adapter 1574666 and extension adapter TITAN 1663911; soft and detachable</p>   | <p><b>1574668</b></p> | <p><b>2</b></p> |
| <p><b>Additional pad II pair</b><br/>                     to attach to TRUMPF coupling point system; used with supporting pads universal; to extend the operating table top; Comfort Plus pad integrally foamed, soft and detachable; Velcro strap; stainless steel load-bearing structure</p>  | <p><b>1594041</b></p> | <p><b>2</b></p> |



| Περιγραφή   | Κωδικός        | Τεμάχια  |
|---|----------------|----------|
| <p><b>Extension adapter</b><br/>adapter for table tops of the TruSystem 7500 and 5500;<br/>for the adaptation of extension accessories; carbon and<br/>stainless steel load-bearing structure</p>  | <b>1574666</b> | <b>2</b> |







| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <p><b>Universal support with joint</b><br/>                     to slip at the extension adapter 1574666 and extension adapter TITAN 1663911; for mounting leg sections, counter traction posts and the MIS-Hip-Aggregat 1574732; with tiltable and removable joints; stainless steel; adjustment: + / - 60°</p>  | 1574667 | 2       |
| <p><b>Double joint strut, pair</b><br/>                     to attach to TRUMPF coupling point system, extension adapter 1574666 and extension adapter TITAN 1663911; stainless steel</p>    | 1574734 | 2       |
| <p><b>Extension strut short</b><br/>                     to attach to double joint strut, pair 1574734 or 454466; for mounting the adapter spindle traction mechanism 1593159; stainless steel load-bearing structure; length: 19 in (485 mm)</p>  | 1597667 | 4       |



| Περιγραφή   | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Adapter spindle traction mechanism</b><br/>                     for mounting the spindle traction mechanism 1574736;<br/>                     stainless steel</p>    | 1593159 | 4       |
| <p><b>Spindle traction mechanism</b><br/>                     with tiltable joint to attach extension shoes;<br/>                     adjustable via ball and socket joint; adjusting of the<br/>                     traction via handwheel; adaptation of the leg length<br/>                     by means of quick adjustment</p>  | 1574736 | 4       |
| <p><b>Holding device for arm and leg</b><br/>                     with clamp; to attach to arm support one hand and pad<br/>                     for leg positioning; stainless steel</p>    | 1317577 | 4       |




| Περιγραφή   | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Clamp basic</b><br/>to attach to side rails; suitable for rods 0.6 in - 0.8 in (16 mm - 20 mm) in diameter; stainless steel</p>    | 1305173 | 4       |
| <p><b>Pad for leg positioning</b><br/>without clamp; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; pad (LxW): 12 in x 5 in (300 mm x 130 mm)</p>  | 1218843 | 4       |

| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <p><b>Leg support</b><br/>without clamp; height adjustable; with leather strap;<br/>stainless steel load-bearing structure</p>    | 4544737 | 4       |
| <p><b>Clamp radial setting</b><br/>to attach to side rails by frontal placement;<br/>for heavy loads; 360° adjustment; suitable for rods<br/>0.5 in - 0.7 in (15 mm - 18 mm) in diameter;<br/>stainless steel</p>  | 1218804 | 4       |


| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <p><b>Rail clamp for extension strut</b><br/>to attach accessories on the extension struts;<br/>stainless steel</p>   | 1218833 | 4       |
| <p><b>Counter traction post 150</b><br/>pad electrically conductive, soft and detachable;<br/>radiolucent rod;<br/>diameter: 7 in (180 mm)<br/>length: 13,8 in (350 mm)</p>  | 1612706 | 2       |

| Περιγραφή  | Κωδικός               | Τεμάχια         |
|--|-----------------------|-----------------|
| <p><b>Counter traction post for femur</b><br/>                     for interventions in lateral position; counter traction post and perineal bow with integrally foamed and soft pad; height adjustable via crank; stainless steel load-bearing structure;<br/>                     height adjustment: 5 in (125 mm)</p>  | <p><b>4544448</b></p> | <p><b>2</b></p> |
| <p><b>Leather extension shoe right</b><br/>                     for adults; four belts for safety support; to attach to spindle traction mechanism</p>   | <p><b>4544540</b></p> | <p><b>2</b></p> |


| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <p><b>Leather extension shoe left</b><br/>                     for adults; four belts for safety support; to attach to spindle traction mechanism</p>                     | 4544541 | 2       |
| <p><b>Leg holder GOEPEL</b><br/>                     without clamp; pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure</p>  | 4544736 | 2       |



| Περιγραφή   | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Clamp radial setting</b><br/>                     to attach to side rails by frontal placement;<br/>                     for heavy loads; 360° adjustment; suitable for rods<br/>                     0.5 in - 0.7 in (15 mm - 18 mm) in diameter;<br/>                     stainless steel</p>  | 1218804 | 2       |
| <p><b>Base support</b><br/>                     stainless steel</p>   | 1218838 | 2       |
| <p><b>Rail clamp for extension strut</b><br/>                     to attach accessories on the extension struts;<br/>                     stainless steel</p>    | 1218833 | 2       |







| Περιγραφή   | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Clamp radial setting</b><br/>to attach to side rails by frontal placement;<br/>for heavy loads; 360° adjustment; suitable for rods<br/>0.5 in - 0.7 in (15 mm - 18 mm) in diameter;<br/>stainless steel</p>  | 1218804 | 2       |



**4.**

| Περιγραφή   | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Extension trolley universal</b><br/>for positioning and transporting OR-adapters including<br/>accessories; with columns (except: SATURN SM) and<br/>during shuttle use as well as MARS 2.01 - 2.06, TITAN<br/>and SATURN Select; stainless steel</p>  | 1504928 | 2       |

| Περιγραφή  | Κωδικός               | Τεμάχια         |
|--|-----------------------|-----------------|
| <p><b>Shoulder positioner</b><br/>           to attach to TRUMPF coupling point system; includes three padded segments; lateral pads are height adjustable and removable; pads integrally foamed, electrically conductive and soft; stainless steel load-bearing structure</p>  | <p><b>4544103</b></p> | <p><b>2</b></p> |
| <p><b>Retainer for head rests</b><br/>           to attach to adapter for head positioning, shoulder segment, shoulder adapter and ophthalmology adapter; stainless steel</p>   | <p><b>1317948</b></p> | <p><b>2</b></p> |


| Περιγραφή  | Κωδικός               | Τεμάχια         |
|--|-----------------------|-----------------|
| <p><b>Connection holder basic asymmetric</b><br/>                     to attach to retainer for head rest; stainless steel</p>    | <p><b>1305449</b></p> | <p><b>2</b></p> |
| <p><b>Head rest for shoulder operation</b><br/>                     radiolucent U-shaped helmet with soft pads and two fixation straps with Velcro fasteners for safe positioning and lateral fixation of the head while patient is anaesthetized; fits various head sizes, ranges from adult to children sizes; removable pads</p>  | <p><b>1296096</b></p> | <p><b>2</b></p> |

| Περιγραφή   | Κωδικός               | Τεμάχια         |
|---|-----------------------|-----------------|
| <p><b>6. Arthroscopy positioning device basic</b><br/>                     without clamp; height adjustable; adjustable to leg width; pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure; height adjustment: 18 in (460 mm)<br/>                     width adjustment: 4 in - 7 in (110 mm - 180 mm)</p>  | <p><b>4544534</b></p> | <p><b>2</b></p> |
| <p><b>Clamp radial setting</b><br/>                     to attach to side rails by frontal placement; for heavy loads; 360° adjustment; suitable for rods 0.5 in - 0.7 in (15 mm - 18 mm) in diameter; stainless steel</p>   | <p><b>1218804</b></p> | <p><b>2</b></p> |


| 7. Περιγραφή   | Κωδικός               | Τεμάχια         |
|--|-----------------------|-----------------|
| <p><b>Counter traction post for tibia</b><br/>                     to attach to the Extension unit or with the Universal seat OT 1576814 to operating table tops; stainless steel</p>  | <p><b>1576815</b></p> | <p><b>2</b></p> |
| <p><b>Adapter counter traction post pivoted</b><br/>                     to attach to the counter traction post for tibia 1576815; stainless steel</p>                               | <p><b>1576816</b></p> | <p><b>2</b></p> |


# ΠΑΠΑΠΟΣΤΟΛΟΥ

ΚΕΝΤΡΟ ΙΑΤΡΙΚΟΥ ΕΞΟΠΛΙΣΜΟΥ

| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <b>Counter traction post 60</b><br>diameter: 2 in (60 mm)<br>length: 10 in (270 mm)<br> | 4544453 | 2       |



8.

| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <b>Arm support table</b><br>with clamp; pad integrally foamed, electrically<br>conductive and soft; stainless steel load-bearing<br>structure;<br>dimensions (LxW): 26 in x 10 in (650 mm x 260 mm)<br> | 4544772 | 2       |


| 9. Περιγραφή  | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Hand operating table</b><br/>with clamp; table top and supports on wheels height adjustable; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; pad (LxW): 32 in x 14 in (820 mm x 360 mm) height adjustment:<br/>table top: 3 in (70 mm)<br/>support on wheels: 16 in (400 mm)<br/>pivot range: +/- 45°</p>  | 4549735 | 2       |


10.

Περιγραφή

| Κωδικός  | Τεμάχια |
|--|---------|
| <p><b>Traction device for humerus</b><br/>without clamp; adjustable to arm length;<br/>stainless steel</p>    | 2       |
| <p><b>Counter traction for humerus</b><br/>without clamp; clamp provides adjustment of height to<br/>match humerus length; padded and radiolucent; stainless<br/>steel load-bearing structure</p>  | 2       |



| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <p><b>WEINBERGER hand traction device</b><br/>                     with clamp; to attach to 4544766 traction device for humerus and 1218832 and 1574736 spindle traction mechanism</p>  | 1218836 | 2       |

| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <p><b>11. Condylen fixation device</b><br/>                     to attach to 1576815 counter traction post for tibia in combination with 1576816 adapter counter traction post pivoted or 1643382 adapter for counter traction post; for exact knee positioning; pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; height adjustment: 5 in (125 mm)</p>  | 1434382 | 2       |

12.

Περιγραφή

**Arthroscopy positioning device basic**

without clamp; height adjustable; adjustable to leg width; pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure; height adjustment: 18 in (460 mm)  
width adjustment: 4 in - 7 in (110 mm - 180 mm)



Κωδικός


4544534



Τεμάχια




1

**13.**

**Περιγραφή**

|  |                                   |                             |
|--|-----------------------------------|-----------------------------|
| <p><b>Spinal cord positioning device</b><br/>to attach to TRUMPF coupling point system at foot end;<br/>used to position patient knees during spinal column<br/>surgeries; with crank for height adjustment; motorized<br/>leg section joints for up and down movements; radiolucent<br/>via C-arm; pad electrically conductive, soft and<br/>detachable; Velcro-fastened pad; stainless steel<br/>load-bearing structure;<br/>height adjustment: 12 in (300 mm)</p>  | <p><b>Κωδικός</b><br/>1361398</p> | <p><b>Τεμάχια</b><br/>2</p> |
|--|-----------------------------------|-----------------------------|

| Περιγραφή   | Κωδικός | Τεμάχια |
|---|---------|---------|
| <p><b>Iliac roll with lateral support</b><br/>                     without clamp; adjustable to table top width; lateral supports can be adjusted to body width; pads integrally foamed, electrically conductive and soft; stainless steel;<br/>                     height adjustment: 16 in (400 mm)</p>   | 4544761 | 2       |
| <p><b>Lateral support height adjustable</b><br/>                     with clamp; to attach to side rails; pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure;<br/>                     pad (LxWxH): 8 in x 4 in x 2 in<br/>                     (215 mm x 100 mm x 50 mm)<br/>                     height adjustment: 2 in (50 mm)<br/>                     pad rotation: 360°</p>  | 1496430 | 4       |

| Περιγραφή  | Κωδικός | Τεμάχια |
|--|---------|---------|
| <p><b>Cushion for sternum support large</b><br/>                     to support the patient's thorax; foam pad;<br/>                     dimensions (LxWxH): 8 in x 23 in x 8 in<br/>                     (200 mm x 580 mm x 200 mm)</p>  | 1436370 | 2       |
| <p><b>Pad for prone position large</b><br/>                     dimensions (LxWxH): 11 in x 9 in x 6 in<br/>                     (280 mm x 230 mm x 145 mm)</p>    | 1237192 | 2       |
| <p><b>Trolley spinal cord positioning device</b><br/>                     for positioning and transport of spinal cord positioning<br/>                     devices; during column or shuttle use; stainless steel</p>                  | 1432336 | 2       |