










8.1 ΣΕΤ ΕΞΑΡΤΗΜΑΤΩΝ ΟΡΘΟΠΑΙΔΙΚΗΣ


Περιγραφή	Κωδικός	Τεμάχια
<p>1.</p> <p>Pad for operation at intervertebral disc foam pad is electrically conductive; dimensions (LxWxH): 33 in x 18 in x 4 in (840 mm x 460 mm x 100 mm)</p> 	1409494	2
<p>2.α</p> <p>Adapter for head positioning to attach to TRUMPF coupling point system; for attaching head positioning accessories; stainless steel; spacer: 9 in (220 mm)</p> 	4544785	2



Περιγραφή	Κωδικός	Τεμάχια
<p>Skull clamp to attach to 1218969 skull clamp adapter; stable and precise three-point skull fixation; cast aluminum</p> 	1218968	2
<p>Basic unit neuro for attaching 1218969 Skull clamp adapter; simultaneous locking of longitudinal and vertical positions; used for prone and supine position of the patient; cast aluminum</p> 	1218970	2

Περιγραφή	Κωδικός	Τεμάχια
<p>Head rest for adults adjustable to head size; to connect to 1218969 skull clamp adapter; with detachable, elastic gel pad; used for prone and supine position of the patient; cast aluminum</p> 	1218974	2
<p>Head rest for children adjustable to head size; to connect to 1218969 skull clamp adapter; with detachable, elastic gel pad; used for prone and supine position of the patient; cast aluminum</p> 	1218973	2



Περιγραφή	Κωδικός	Τεμάχια
<p>Skull clamp adapter connection between 1218970 basic unit neuro, 1218968 skull clamp or 4145424 multi-purpose head rest; cast aluminum</p> 	1218969	2
<p>Skull pins for adults sterilisable; made of stainless steel; package with three pieces</p> 	4145420	2


Περιγραφή	Κωδικός	Τεμάχια
Skull pins for children sterilisable; made of stainless steel; package with three pieces 	4145421	2



2β. Retainer for head rests to attach to adapter for head positioning, shoulder segment, shoulder adapter and ophthalmology adapter; stainless steel 	Κωδικός 1317948	Τεμάχια 2
---	--------------------	--------------




Περιγραφή	Κωδικός	Τεμάχια
<p>Connection holder basic asymmetric to attach to retainer for head rest; stainless steel</p> 	1305449	2
<p>Head rest basic with clamp; to attach to connection holder; Comfort Plus pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure; dimensions (LxV): 10 in x 8 in (250 mm x 200 mm)</p> 	4549730	2



3.



Περιγραφή	Κωδικός	Τεμάχια
<p>Supporting pad universal pad for extension adapter 1574666 and extension adapter TITAN 1663911; soft and detachable</p> 	<p>1574668</p>	<p>2</p>
<p>Additional pad II pair to attach to TRUMPF coupling point system; used with supporting pads universal; to extend the operating table top; Comfort Plus pad integrally foamed, soft and detachable; Velcro strap; stainless steel load-bearing structure</p> 	<p>1594041</p>	<p>2</p>



Περιγραφή	Κωδικός	Τεμάχια
<p>Extension adapter adapter for table tops of the TruSystem 7500 and 5500; for the adaptation of extension accessories; carbon and stainless steel load-bearing structure</p> 	1574666	2



Περιγραφή	Κωδικός	Τεμάχια
<p>Universal support with joint to slip at the extension adapter 1574666 and extension adapter TITAN 1663911; for mounting leg sections, counter traction posts and the MIS-Hip-Aggregat 1574732; with tiltable and removable joints; stainless steel; adjustment: + / - 60°</p> 	1574667	2
<p>Double joint strut, pair to attach to TRUMPF coupling point system, extension adapter 1574666 and extension adapter TITAN 1663911; stainless steel</p> 	1574734	2
<p>Extension strut short to attach to double joint strut, pair 1574734 or 454466; for mounting the adapter spindle traction mechanism 1593159; stainless steel load-bearing structure; length: 19 in (485 mm)</p>	1597667	4



Περιγραφή	Κωδικός	Τεμάχια
<p>Adapter spindle traction mechanism for mounting the spindle traction mechanism 1574736; stainless steel</p> 	1593159	4
<p>Spindle traction mechanism with tiltable joint to attach extension shoes; adjustable via ball and socket joint; adjusting of the traction via handwheel; adaptation of the leg length by means of quick adjustment</p> 	1574736	4
<p>Holding device for arm and leg with clamp; to attach to arm support one hand and pad for leg positioning; stainless steel</p> 	1317577	4




Περιγραφή	Κωδικός	Τεμάχια
<p>Clamp basic to attach to side rails; suitable for rods 0.6 in - 0.8 in (16 mm - 20 mm) in diameter; stainless steel</p> 	1305173	4
<p>Pad for leg positioning without clamp; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; pad (LxW): 12 in x 5 in (300 mm x 130 mm)</p> 	1218843	4


Περιγραφή	Κωδικός	Τεμάχια
<p>Leg support without clamp; height adjustable; with leather strap; stainless steel load-bearing structure</p> 	4544737	4
<p>Clamp radial setting to attach to side rails by frontal placement; for heavy loads; 360° adjustment; suitable for rods 0.5 in - 0.7 in (15 mm - 18 mm) in diameter; stainless steel</p> 	1218804	4

Περιγραφή	Κωδικός	Τεμάχια
<p>Rail clamp for extension strut to attach accessories on the extension struts; stainless steel</p> 	1218833	4
<p>Counter traction post 150 pad electrically conductive, soft and detachable; radiolucent rod; diameter: 7 in (180 mm) length: 13,8 in (350 mm)</p> 	1612706	2


Περιγραφή	Κωδικός	Τεμάχια
<p>Counter traction post for femur for interventions in lateral position; counter traction post and perineal bow with integrally foamed and soft pad; height adjustable via crank; stainless steel load-bearing structure; height adjustment: 5 in (125 mm)</p> 	4544448	2
<p>Leather extension shoe right for adults; four belts for safety support; to attach to spindle traction mechanism</p> 	4544540	2



Περιγραφή	Κωδικός	Τεμάχια
<p>Leather extension shoe left for adults; four belts for safety support; to attach to spindle traction mechanism</p> 	4544541	2
<p>Leg holder GOEPEL without clamp; pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure</p> 	4544736	2



Περιγραφή	Κωδικός	Τεμάχια
<p>Clamp radial setting to attach to side rails by frontal placement; for heavy loads; 360° adjustment; suitable for rods 0.5 in - 0.7 in (15 mm - 18 mm) in diameter; stainless steel</p> 	1218804	2
<p>Base support stainless steel</p> 	1218838	2
<p>Rail clamp for extension strut to attach accessories on the extension struts; stainless steel</p> 	1218833	2



Περιγραφή	Κωδικός	Τεμάχια
<p>Clamp radial setting to attach to side rails by frontal placement; for heavy loads; 360° adjustment; suitable for rods 0.5 in - 0.7 in (15 mm - 18 mm) in diameter; stainless steel</p> 	1218804	2



4. Περιγραφή


Περιγραφή	Κωδικός	Τεμάχια
<p>Extension trolley universal for positioning and transporting OR-adapters including accessories; with columns (except: SATURN SM) and during shuttle use as well as MARS 2.01 - 2.06, TITAN and SATURN Select; stainless steel</p> 	1504928	2

Περιγραφή	Κωδικός	Τεμάχια
<p>Shoulder positioner to attach to TRUMPF coupling point system; includes three padded segments; lateral pads are height adjustable and removable; pads integrally foamed, electrically conductive and soft; stainless steel load-bearing structure</p> 	<p>4544103</p>	<p>2</p>
<p>Retainer for head rests to attach to adapter for head positioning, shoulder segment, shoulder adapter and ophthalmology adapter; stainless steel</p> 	<p>1317948</p>	<p>2</p>


Περιγραφή	Κωδικός	Τεμάχια
<p>Connection holder basic asymmetric to attach to retainer for head rest; stainless steel</p> 	1305449	2
<p>Head rest for shoulder operation radiolucent U-shaped helmet with soft pads and two fixation straps with Velcro fasteners for safe positioning and lateral fixation of the head while patient is anaesthetized; fits various head sizes, ranges from adult to children sizes; removable pads</p> 	1296096	2


Περιγραφή	Κωδικός	Τεμάχια
<p>6. Arthroscopy positioning device basic without clamp; height adjustable; adjustable to leg width; pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure; height adjustment: 18 in (460 mm) width adjustment: 4 in - 7 in (110 mm - 180 mm)</p> 	<p>4544534</p>	<p>2</p>
<p>Clamp radial setting to attach to side rails by frontal placement; for heavy loads; 360° adjustment; suitable for rods 0.5 in - 0.7 in (15 mm - 18 mm) in diameter; stainless steel</p> 	<p>1218804</p>	<p>2</p>

Περιγραφή	Κωδικός	Τεμάχια
<p>Counter traction post for tibia to attach to the Extension unit or with the Universal seat OT 1576814 to operating table tops; stainless steel</p> 	1576815	2
<p>Adapter counter traction post pivoted to attach to the counter traction post for tibia 1576815; stainless steel</p> 	1576816	2

Περιγραφή	Κωδικός	Τεμάχια
Counter traction post 60 diameter: 2 in (60 mm) length: 10 in (270 mm) 	4544453	2



8.


Περιγραφή	Κωδικός	Τεμάχια
Arm support table with clamp; pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure; dimensions (LxW): 26 in x 10 in (650 mm x 260 mm) 	4544772	2


9. Περιγραφή	Κωδικός	Τεμάχια
<p>Hand operating table with clamp; table top and supports on wheels height adjustable; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; pad (LxW): 32 in x 14 in (820 mm x 360 mm) height adjustment: table top: 3 in (70 mm) support on wheels: 16 in (400 mm) pivot range: +/- 45°</p> 	4549735	2

10.

Περιγραφή

Περιγραφή	Κωδικός	Τεμάχια
<p>Traction device for humerus without clamp; adjustable to arm length; stainless steel</p> 	<p>4544766</p>	<p>2</p>
<p>Counter traction for humerus without clamp; clamp provides adjustment of height to match humerus length; padded and radiolucent; stainless steel load-bearing structure</p> 	<p>4544765</p>	<p>2</p>

Περιγραφή	Κωδικός	Τεμάχια
<p>WEINBERGER hand traction device with clamp; to attach to 4544766 traction device for humerus and 1218832 and 1574736 spindle traction mechanism</p> 	1218836	2

Περιγραφή	Κωδικός	Τεμάχια
<p>11. Condylen fixation device to attach to 1576815 counter traction post for tibia in combination with 1576816 adapter counter traction post pivoted or 1643382 adapter for counter traction post; for exact knee positioning; pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; height adjustment: 5 in (125 mm)</p> 	1434382	2

12.

Περιγραφή

Arthroscopy positioning device basic

without clamp; height adjustable; adjustable to leg width; pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure; height adjustment: 18 in (460 mm)
width adjustment: 4 in - 7 in (110 mm - 180 mm)



Κωδικός

4544534

Τεμάχια

1

13.

Περιγραφή

Spinal cord positioning device

to attach to TRUMPF coupling point system at foot end; used to position patient knees during spinal column surgeries; with crank for height adjustment; motorized leg section joints for up and down movements; radiolucent via C-arm; pad electrically conductive, soft and detachable; Velcro-fastened pad; stainless steel load-bearing structure; height adjustment: 12 in (300 mm)








Κωδικός

1361398

Τεμάχια

2

Περιγραφή	Κωδικός	Τεμάχια
<p>Iliac roll with lateral support without clamp; adjustable to table top width; lateral supports can be adjusted to body width; pads integrally foamed, electrically conductive and soft; stainless steel; height adjustment: 16 in (400 mm)</p> 	4544761	2
<p>Lateral support height adjustable with clamp; to attach to side rails; pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; pad (LxWxH): 8 in x 4 in x 2 in (215 mm x 100 mm x 50 mm) height adjustment: 2 in (50 mm) pad rotation: 360°</p> 	1496430	4

Περιγραφή	Κωδικός	Τεμάχια
<p>Cushion for sternum support large to support the patient's thorax; foam pad; dimensions (LxWxH): 8 in x 23 in x 8 in (200 mm x 580 mm x 200 mm)</p> 	1436370	2
<p>Pad for prone position large dimensions (LxWxH): 11 in x 9 in x 6 in (280 mm x 230 mm x 145 mm)</p> 	1237192	2
<p>Trolley spinal cord positioning device for positioning and transport of spinal cord positioning devices; during column or shuttle use; stainless steel</p> 	1432336	2