






8.4 – 8.5 ΣΕΤ ΕΞΑΡΤΗΜΑΤΩΝ ΓΥΝΑΙΚΟΛΟΓΙΚΗΣ – ΟΥΡΟΛΟΓΙΚΗΣ


1.



Περιγραφή	Κωδικός	Τεμάχια
Drainage bowl with sieve; suitable for all drainage bowl adapters; hose connectors for drainage; stainless steel 	4544618	2
Drainage bowl adapter with clamp; to attach to side rail; guide rail allows precise adjustment in x- and y-directions; stainless steel 	4544619	2

2.



Περιγραφή	Κωδικός	Τεμάχια
<p>Leg holder GOEPEL without clamp; pad integrally foamed, electrically conductive and soft; stainless steel load-bearing structure</p> 	4544736	4
<p>Clamp radial setting to attach to side rails by frontal placement; for heavy loads; 360° adjustment; suitable for rods 0.5 in - 0.7 in (15 mm - 18 mm) in diameter; stainless steel</p> 	1218804	4

Περιγραφή	Κωδικός	Τεμάχια
3. Kick bucket without pushing handle with stainless steel bucket; dimensions (LxWxH): 18 in x 18 in x 15 in (470 mm x 470 mm x 380 mm) 	1218658	2

Περιγραφή	Κωδικός	Τεμάχια
4. Lateral support 215 to attach to lateral supports; pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; pad (LxWxH): 4 in x 8 in x 2 in (100 mm x 215 mm x 50 mm) 	4544721	2

Περιγραφή	Κωδικός	Τεμάχια
<p>Universal lateral support height adjusta with clamp; to attach to side rails by frontal placement; to attach lateral supports; pad can be individually adjusted; stainless steel; height adjustment: 6 in (150 mm)</p> 	<p>4544719</p>	<p>2</p>
<p>Pelvis support with clamp; for lateral positioning; pads can be adjusted to pelvis width; pads integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; pad (LxWxH): 3 in x 3 in x 2 in (85 mm x 85 mm x 50 mm)</p> 	<p>4544773</p>	<p>2</p>

5.

Περιγραφή	Κωδικός	Τεμάχια
<p>Leg holder left without clamp; easily adjustable through gas-spring supported mechanism; one handed operating; marks for position adjustment; only with leg holder right 1228197 available (as a pair)</p> 	1228196	2
<p>Leg holder right without clamp; easily adjustable through gas-spring supported mechanism; one handed operating; marks for position adjustment; only with leg holder left 1228196 available (as a pair)</p> 	1228197	2

ΠΑΠΑΠΟΣΤΟΛΟΥ

ΚΕΝΤΡΟ ΙΑΤΡΙΚΟΥ ΕΞΟΠΛΙΣΜΟΥ

Περιγραφή	Κωδικός	Τεμάχια
Clamp for leg holder for legholder left 1228196 and leg holder right 1228197 	1228198	4