








1<sup>η</sup> Ομάδα: Χειρουργικές Τράπεζες (τεμάχια 21)




Κάθε μία από τις 21 χειρουργικές τράπεζες, περιλαμβάνει τα κάτωθι:

Περιγραφή	Κωδικός	Τεμάχια
<p><b>Mobile column TruSystem 7500</b> mobile operating table column with electrically motorized adjustment of height, lateral tilting and Trendelenburg/anti-Trendelenburg; power supply integrated in column as well as integrated high-capacity battery; column and table top can be operated with the column keypad or the cable or wireless remote control or foot control; column cladding and base plate are made of stainless steel; with NanoProtect</p> 	1717023	1

Περιγραφή	Κωδικός	Τεμάχια
<p><b>Operating table top U24 W</b> washable operating table top; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable</p> 	1572671	2
<p><b>Head section COMFORT</b> to attach to TRUMPF coupling point system; with single joint and one hand operation; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; with NanoProtect; dimensions (LxW): 11 in x 20 in (280 mm x 500 mm) pad (LxW): 9 in x 15 in (225 mm x 375 mm) adjustment range: +/- 55°</p> 	1511441	2

Περιγραφή	Κωδικός	Τεμάχια
<p><b>Upper back section</b>                      to attach to TRUMPF coupling point system; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; with NanoProtect; length: 9 in (225 mm)</p> 	1424811	2
<p><b>Leg section two parts spread joint</b>                      to attach to TRUMPF coupling point system; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable; stainless steel load-bearing structure; with NanoProtect; length: 29 in (740 mm)</p> 	1481482	2

Περιγραφή	Κωδικός	Τεμάχια
<p><b>Cassette rails upper back section</b>                      installation-kit for upper back section and pelvis extension SATURN</p> 	4544143	2
<p><b>Shuttle 3.6 360</b>                      for transporting just the operating table top or the table top plus a mobile operating table column; with continuous Trendelenburg/ anti-Trendelenburg adjustment via crank and height adjustment via pedal; running gear with four electrically conductive swivel wheels, pedal for braking and a 5th wheel to support directional travel; made of aluminium;                      Trendelenburg/ anti-Trendelenburg: +/- 15°                      height adjustment: 9 in (230 mm)                      maximum load: 794 lbs (360 kg)</p> 	1489491	2


Περιγραφή	Κωδικός	Τεμάχια
<b>Cable remote control TS7500</b> Attachable 	1633968	1
<b>Remote control wireless TS 7500</b> Rechargeable 	1633967	1
<b>Charging unit mobile AC TS7500</b> for operating table system TruSystem 7500 and 5500; 100-240 VAC 	1557243	1

# ΠΑΠΑΠΟΣΤΟΛΟΥ

ΚΕΝΤΡΟ ΙΑΤΡΙΚΟΥ ΕΞΟΠΛΙΣΜΟΥ

Ένας από τους κάτωθι ποδοδιακόπτες, κατ' επιλογή σας

Περιγραφή	Κωδικός	Τεμάχια
<b>Foot switch HKN</b> to control height, Trendelenburg/ anti-Trendelenburg and tilt; case made of aluminium	1258883	
<b>Foot switch HBN</b> to control height, Trendelenburg/ anti-Trendelenburg and leg section; case made of aluminium	1258884	
<b>Foot switch HRN</b> to control height, Trendelenburg/ anti-Trendelenburg and back section; case made of aluminium	1258885	
<b>Foot switch HLN</b> to control height, Trendelenburg/ anti-Trendelenburg and longitudinal shift; case made of aluminium	1258887	
<b>Foot switch LBR</b> to control longitudinal shift, leg section and back section; case made of aluminium	1258888	1
<b>Foot switch HBK</b> to control height, leg section and tilt; case made of aluminium	1364476	
<b>Foot switch HRK</b> to control height, back section and tilt; case made of aluminium	1364477	





# ΠΑΠΑΠΟΣΤΟΛΟΥ



ΚΕΝΤΡΟ ΙΑΤΡΙΚΟΥ ΕΞΟΠΛΙΣΜΟΥ

Περιγραφή


Κωδικός

Τεμάχια

<p><b>Arm support 450</b> with clamp; to attach to side rails by frontal placement; swivel mounted and height adjustable; Comfort Plus pad integrally foamed, electrically conductive, soft and detachable; black polyester Velcro strap; stainless steel load-bearing structure; dimensions (LxW): 24 in x 6 in (600 mm x 140 mm) height adjustment: 7 in (170 mm)</p> 	<p>4549778</p>	<p>2</p>
<p><b>Body strap with locking clamps</b> to attach to side rail; Velcro fastening system provides wide range of adjustment for abdominal girth; resistant to disinfectants and washable; black polyester strap; stainless steel locking clamp; dimensions (LxW): 71 in x 4 in (1800 mm x 100 mm)</p> 	<p>4544702</p>	<p>1</p>

Περιγραφή	Κωδικός	Τεμάχια
<p><b>Infusion pole with 4 hooks</b> with clamp; to attach to side rails by frontal placement; with ball joint; height adjustable; stainless steel; height adjustment: 24 in (620 mm)</p> 	4544708	2
<p><b>Anaesthesia screen</b> without clamp; to attach to side rails; electrically insulated plus a safety ring; height adjustable via clamp; stainless steel; vertical rod diameter: 0.6 in (16 mm) height adjustment: 20 in (500 mm)</p> 	1218820	1



Περιγραφή	Κωδικός	Τεμάχια
<p><b>Y-clamp</b> to attach to side rail by frontal placement; suitable for rods 0.5 in - 0.7 in (14 mm - 18 mm) in diameter; stainless steel</p> 	<b>4544701</b>	<b>1</b>