








3.3 Στήλη Χειρουργού (τεμάχια 6)

Κάθε μία στήλη χειρουργού από τις έξι προσφερόμενες περιλαμβάνει τα κάτωθι:

Περιγραφή	Κωδικός	Τεμάχια
System flange tandem ML/HL2000 LK560 35,0 kg 	0337323	1

Περιγραφή	Κωδικός	Τεμάχια
<p>Ceiling fixation set H01/8 ZD 1200 for intermediate ceiling max. distance intermediate/bare ceiling 1200 mm consisting of: 1x Ceiling anchor plate 560/605/ M16/M22 8x Distance unit M16/M22</p> 	<p>0337447</p>	<p>1</p>
<p>Interface plate standard with interfaces</p> 	<p>0337311</p>	<p>1</p>

Περιγραφή	Κωδικός	Τεμάχια
<p>Ceiling cover Tandem 640 x 640 x 100</p> 	0337174	1
<p>Cover blank ceiling cover / panel Tandem</p> <p>PHOTO NOT AVAILABLE</p>	0337039	1
<p>TruPort 7500 - 2300</p> <p>2 swivel arms 1000 + 1300 mm motor-driven vertical movement lift speed 30 mm/s</p> 	1638923	1

Περιγραφή	Κωδικός	Τεμάχια
<p>TruPort M system head System head dimensions 310 mm x 220 mm</p> 	1636193	1
<p>System head length 1000</p> 	1636255	1


ΠΑΠΑΠΟΣΤΟΛΟΥ



ΚΕΝΤΡΟ ΙΑΤΡΙΚΟΥ ΕΞΟΠΛΙΣΜΟΥ


Περιγραφή


Κωδικός

Τεμάχια

<p>Elektropneum. brake HL-top, middle joint suitable for TruPort 7000 units PHOTO NOT AVAILABLE</p>	<p>1636190</p>	<p>1</p>
<p>Elektropneum.brake in head column joint suitable for TruPort 5000 / TruPort 7000 not suitable for TruPort x500 not suitable for TruPort x510 PHOTO NOT AVAILABLE</p>	<p>1636191</p>	<p>1</p>
<p>Shelf 530x480 ComfortLine incl. ComfortLine interface TruPort shelf 530 mm wide, 480 mm depth</p> 	<p>1636318</p>	<p>2</p>

Περιγραφή	Κωδικός	Τεμάχια
<p>Operating handle shelf Additional costs for shelf with ComfortLine interface</p> 	1636351	1
<p>Double swivel arm 240 mm - VESA 75/100 Monitor holder incl. Adapter plate VESA 75/100 Lengths of swivel arms 240 mm + 240 mm Fixation with - Adapter lateral - Fixation clamp</p> <p>2,5 kg</p> 	0337432	1

Περιγραφή	Κωδικός	Τεμάχια
<p>Adapter monitor holder ComfortLine incl. ComfortLine interface suitable for all monitor holders</p> <p>1,1 kg</p> 	1640193	1
<p>Circuit cable/wiring / PowerLine 1 x circuit cable 3 x 2,5qmm max. 20A per support head 4 circuit cable required, from support head length 1000 also 8 circuit cable possible</p> <p>PHOTO NOT AVAILABLE</p>	1638285	4
<p>Power socket / PowerLine incl. wiring inside supply modul</p> <p>PHOTO NOT AVAILABLE</p>	1638283	12

Περιγραφή	Κωδικός	Τεμάχια
<p>Earthing pin in following versions available: module S: 4 items module L: 8 items</p> 	1638286	12
<p>Gas outlet medical air 3-5 (free issue) Installation PHOTO NOT AVAILABLE</p>	0337180	2
<p>Gas outlet medical air 7-10 (free issue) Installation PHOTO NOT AVAILABLE</p>	0337185	2
<p>Gas outlet medical air 7-10 (free issue) for <u>neurosurgery operating rooms</u> Installation PHOTO NOT AVAILABLE</p>	0337185	2

Περιγραφή	Κωδικός	Τεμάχια
Gas outlet oxygen (free issue) Installation PHOTO NOT AVAILABLE	0337179	1
Gas outlet vacuum (free issue) Installation PHOTO NOT AVAILABLE	0337182	2
Low voltage cut out for preperation acc. the documentation cost for each cut out PHOTO NOT AVAILABLE	1714358	4
Cover Low voltage cut out Plastic or aluminium costs for each cover PHOTO NOT AVAILABLE	1714359	4
Cable conduit AD 28,5 (ID 20,5) with pull wire from system flange to system head PHOTO NOT AVAILABLE	0337099	1